

SCALP VEIN SET

CHIRAFLEX

Instruction for Use Scalp vein set for Single Use



CONTAINS or PRESENCE of PHTHALATES: Bis (2-ethylhexyl)phthalate (DEHP) 

- Description:** Scalp vein set set for single use. Sterile, Non-toxic, Non-pyrogenic.
- Components:** The products composed of cannula, handle, soft tube, hub and protector.
- Application scope:** This product is used in routine blood collection procedures and/or the short-term infusion of intravenous fluids. The winged needle is designed with a safety shield, which can be activated to cover the needle immediately following blood collection/infusion to aid in the protection against accidental needlestick injury. The product is to be used by appropriately trained healthcare professionals only in accordance with these instructions.
- Prohibition:** ---



1. Before usage, confirm the sterility validity, packing completeness, the lot No. and expiry date are in order;
2. The products should be used immediately after envelope opened, discard after use;
3. It is not allowed for use if envelope damaged, cap falling off and validity expired;
4. If there is unusual cased during transfusion/infusion, treat according to the clinical operation procedures;
5. This device contains phthalates DEHP;
6. The infusion of lipid based solution could increase the risk of endocrin disruption in children, infants or pregnant women should be avoided;
7. Re-sterilization or re-use could cause infection or cross-interactions.



1. Before usage, open the envelope from the „open here“ labeled place;
2. Examine individual package for integrity of packaging prior to use. If packaging has been torn, do not use;
3. Any used needle is considered contaminated and should be disposed of in an approved “sharps“ biohazard container immediately after use;
4. Keep hands behind needle at all times during use and disposal.

| | | | |
|---|------------------------------|--|---------------------------------------|
|  | Attention |  | Manufacture by |
|  | See instruction for use |  | LATEX Free product |
|  | Product reference/Art. No. |  | Do not use if unit package is damaged |
|  | For single use only |  | Non pyrogenic |
|  | Batch Number |  | Keep away from sunlight |
|  | Date of Manufacturing |  | Temperature Limitation |
|  | Use by / Expiry date |  | Keep away from Rain |
|  | Sterilised by Ethylene Oxide | | |


CHIRANA T. Injecta, a.s.
 Nám. Dr. Schweitzera 194
 916 01 Stará Turá
 Slovak Republic
 Tel.: +421/32/7709 991
 Fax.: +421/32/7709 956
www.t-injecta.com