

CHIRANA® ENTERAL syringes



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ 60 ml enteral syringes with female ENFit connector (ISO 80369-3)
- ✓ ENFit connector provides a simple way to reduce the risk of enteral tube feeding misconnections and improve patient safety
- ✓ Does not allow connectivity with any other connector for any other clinical use
- ✓ Provides a locking feature that signals the appropriate connection and stays in place
- ✓ Has a female connector end for any administration sets and syringes that fit into a male patient-access feeding tube port
- ✓ Made of polypropylene and lubricated synthetic rubber

Size	Pieces per Box	Pieces per Case	Item Code
60 ml Enteral syringe	50	200	CH03050E