

CHIRANA® 3-part PERFUSION syringes



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ manufactured and intended for PERFUSION applications
- ✓ compatible with perfusion pumps of world leader manufacturers
- ✓ compatible with dialysis units of world leader manufacturers
- ✓ black graduation, easy to read, DOUBLE graduation on barrel and piston
- ✓ made of polypropylene and lubricated synthetic rubber
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ sterile, individually blister packed

Size	Pieces per Case	Item Code
20 ml transparent for dialysis machine use	100	CH03020PT
30 ml transparent for dialysis machine use	100	CH03030PT
50 ml transparent	100	CH03050LL
50 ml transparent	100	CH03050PT

CHIRANA® 3-part PERFUSION syringes

Size	Pieces per Case	Item Code
50 ml opaque	100	CH03050PO
50 ml transparent with needle 2,0x30 mm	100	CH030502030PT
50 ml transparent BB	100	CH03050PTB
50 ml opaque BB	100	CH03050POB
50 ml transparent BB with needle 2,0x30 mm	100	CH030502030PTB
50 ml opaque BB with needle 2,0x30 mm	100	CH030502030POB
50 ml black BB with needle 2,0x30 mm	100	CH030502030PBB
50 ml transparent F	100	CH03050PTF
50 ml opaque F	100	CH03050POF
50 ml transparent F with needle 1,8x40 mm	100	CH030501840PTF

✓ TRANSPARENT, OPAQUE and BLACK barrel material versions are available