

CHIRAVAC™ Blood Collection needles



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ screw-on needle hub compatible with all available holders
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ easy-to-open packaging and manipulation
- ✓ comfortable collection experience
- ✓ sterile, individually packed

| Size metric | Size G/'' | Color code | Pieces per Box | Pieces per Case | Item Code |
|--------------|---------------|------------|----------------|-----------------|------------|
| 0,60 x 38 mm | 23G x 1 1/2'' | blue | 100 | 1000 | CHBCN23112 |
| 0,70 x 38 mm | 22G x 1 1/2'' | black | 100 | 1000 | CHBCN22112 |
| 0,80 x 38 mm | 21G x 1 1/2'' | green | 100 | 1000 | CHBCN21112 |
| 0,90 x 38 mm | 20G x 1 1/2'' | yellow | 100 | 1000 | CHBCN20112 |
| 1,20 x 38 mm | 18G x 1 1/2'' | pink | 100 | 1000 | CHBCN18112 |
| 0,60 x 25 mm | 23G x 1'' | blue | 100 | 1000 | CHBCN23100 |

CHIRAVAC™ Blood Collection needles

| Size metric | Size G/” | Color code | Pieces per Box | Pieces per Case | Item Code |
|--------------|----------|------------|----------------|-----------------|------------|
| 0,70 x 25 mm | 22G x 1” | black | 100 | 1000 | CHBCN22100 |
| 0,80 x 25 mm | 21G x 1” | green | 100 | 1000 | CHBCN21100 |
| 0,90 x 25 mm | 20G x 1” | yellow | 100 | 1000 | CHBCN20100 |
| 1,20 x 25 mm | 18G x 1” | pink | 100 | 1000 | CHBCN18100 |

- ✓ intended for blood collection with vacuum tubes
- ✓ made of polypropylene hub and stainless steel tube
- ✓ thin wall tube material
- ✓ color coded needle hub