

CHIRANA® 3-part disposable syringes with needles



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ ergonomic design, smooth flow, highly transparent barrel, tight
- ✓ black extended graduation, easy to read
- ✓ designed tip for best connection, scratch-free for glove manipulation
- ✓ made of polypropylene and lubricated synthetic rubber
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ sterile, individually blister packed
- ✓ no dead space volume for 1ml
- ✓ Luer Slip tip
- ✓ syringe with needle available in other needle combinations
- ✓ needle is side-packed in single blister package

Size	Pieces per Box	Pieces per Case	Item Code LUER
2 ml with needle 0,60 x 30 mm (23 G x 1 1/4")	100	3000	CH03002L0630
5 ml with needle 0,70 x 30 mm (22 G x 1 1/4")	100	2400	CH03005L0730
10 ml with needle 0,80 x 40 mm (21 G x 1 1/2")	100	1600	CH03010L0840
20 ml with needle 0,80 x 40 mm (21 G x 1 1/2")	50	800	CH03020L0840
50 ml Luer with needle 1,20 x 40 mm (18 G x 1 1/2")	50	300	CH03050L1240