


MEDICAL DEVICES



NEW



Company profile

CTI group is a Constellation of 3 European manufacturers of world-class quality products for the global healthcare sector. It multiplies the strength of CHIRANA T. Injecta, VITREX, and CHIRMAX companies on the global healthcare market.

CTI group is Committed To Innovation and reflects this status to daily team work and serving to customers.

CHIRANA T. Injecta, a.s. as a world producer of first-class medical products and services carries on the tradition of the medical devices production in Stará Turá. The manufacture beginnings go back to 1935 when production of glass syringes was launched. Thanks to the dynamic production programme and continuous development of new products the company has started to manufacture dispos-

able plastic syringes in 1970 followed by disposable needles in 1975. One of the most important company's achievement is dated in 1991 when the production of lancets has been established. Present product range is focused primary in the field of hypodermic syringes, needles, canulas, lancets and other medical products.

The company has implemented Quality Management System in compliance with the requirements of EN ISO 13485:2012 for the design, manufacturing, assembly, sterilization and sales. Products are accredited with CE mark in compliance with the EU Medical Device Directive 93/42/EEC.

We Supply Values for Life

Chirana® MEDOJECT®  CHIRMAX CHIRAVAC™ MEDOJECT® fine





CHIRANA® 2-part disposable syringes

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ ergonomic design, smooth flow, highly transparent barrel, tight
- ✓ black extended graduation, easy to read, double stop lines, guaranteed 0 ml position
- ✓ designed tip for best connection, double stop lines, scratch-free for glove manipulation
- ✓ available with Luer tip
- ✓ made of polypropylene/polyethylene, rubber-free
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ sterile, individually blister packed

Size	Pieces per Box	Pieces per Case	Item Code LUER
2 ml	100	2800	CH002L
5 ml	100	1900	CH005L
10 ml	100	1200	CH010L
20 ml	80	960	CH020L



CHIRANA® 2-part disposable syringes with needles

Size	Pieces per Box	Pieces per Case	Item Code LUER
2 ml with needle 0,60 x 30 mm (23 G x 1 1/4")	100	2800	CH002L0630
5 ml with needle 0,70 x 30 mm (22 G x 1 1/4")	100	1900	CH005L0730
10 ml with needle 0,80 x 40 mm (21 G x 1 1/2")	100	1200	CH010L0840
20 ml with needle 0,80 x 40 mm (21 G x 1 1/2")	80	960	CH020L0840

- ✓ syringe with needle available in other needle combinations
- ✓ needle is side-packed in single blister package



CHIRANA® 3-part disposable syringes

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ ergonomic design, smooth flow, highly transparent barrel, tight
- ✓ black extended graduation, easy to read
- ✓ designed tip for best connection, scratch-free for glove manipulation
- ✓ available with Luer and Luer lock
- ✓ made of polypropylene, polyethylene and lubricated synthetic rubber
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ sterile, individually blister packed

Size	Pieces per Box	Pieces per Case	Item Code LUER	Item Code LUER LOCK
1 ml	100	2000	CH03001L	CH03001LL
2 ml	100	3000	CH03002L	CH03002LL
3 ml	100	3000	CH03003L	CH03003LL
5 ml	100	2400	CH03005L	CH03005LL
10 ml	100	1600	CH03010L	CH03010LL
20 ml	50	800	CH03020L	CH03020LL
30 ml	40	640	CH03030L	CH03030LL
50 ml	50	450	CH03050L	CH03050LL

- ✓ low dead space volume for 1 ml

- ✓ centrally positioned Luer lock tip



CHIRANA® 3-part disposable syringes with needles

Size	Pieces per Box	Pieces per Case	Item Code LUER
2 ml with needle 0,60 x 30 mm (23 G x 1 1/4")	100	3000	CH03002L0630
5 ml with needle 0,70 x 35 mm (22 G x 1 2/5")	100	2400	CH03005L0735
10 ml with needle 0,80 x 40 mm (21 G x 1 1/2")	100	1600	CH03010L0840
20 ml with needle 0,80 x 40 mm (21 G x 1 1/2")	50	800	CH03020L0840
50 ml Luer with needle 1,20 x 40 (18 G x 1 1/2")	50	300	CH03050L1240

- ✓ syringe with needle available in other needle combinations

- ✓ needle is side-packed in single blister package



CHIRANA® PERFUSION 3-part syringes

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ manufactured and intended for PERFUSION applications
- ✓ compatible with perfusion pumps of world leader manufacturers
- ✓ compatible with dialysis units of world leader manufacturers
- ✓ black graduation, easy to read, DOUBLE graduation on barrel and piston
- ✓ made of polypropylene, polyethylene and lubricated synthetic rubber
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ sterile, individually blister packed

Size	Pieces per Case	Item Code
20 ml transparent for dialysis machine use	100	CH03020PT
30 ml transparent for dialysis machine use	100	CH03030PT
50 ml transparent	100	CH03050PT
50 ml opaque	100	CH03050PO
50 ml transparent with needle 2,0x30 mm	100	CH030502030PT
50 ml transparent BB	100	CH03050PTB
50 ml opaque BB	100	CH03050POB
50 ml transparent BB with needle 2,0x30 mm	100	CH030502030PTB
50 ml transparent F	100	CH03050PTF
50 ml opaque F	100	CH03050POF
50 ml transparent F with needle 1,8x40 mm	100	CH030501840PTF

- ✓ TRANSPARENT, OPAQUE and BLACK barrel material versions are available
- ✓ available with aspiration needle and filter needle
- ✓ possibility of ETO and GAMMA sterilisation

CHIRANA® 3-part disposable Catheter syringes

Size	Pieces per Box	Pieces per Case	Item Code LUER
50 ml Catheter		100	CH03050C
50 ml Catheter with Luer adapter		100	CH03050CLA
100 ml Catheter with Luer adapter	25	100	CH03100C
150 ml Catheter Janette®		50	CH03150CS

- ✓ intended for irrigation application
- ✓ 150 ml Catheter Janette® under declared conditions suitable for autoclavable process



CHIRANA® Insulin syringes

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ needle permanently fixed to syringe body
- ✓ no dead space volume, Insulin and Tuberculin graduation
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ sterile, individually blister packed

Type	Size	Size G/"	Size metric	Pieces per Box	Pieces per Case	Item Code
INSULIN U40	1 ml	29G x 1/2"	12,7 mm	100	2000	CHINS4129
		27G x 1/2"	12,7 mm	100	2000	CHINS0127
		29G x 5/16"	8 mm	100	2000	CHINS01298
		29G x 1/2"	12,7 mm	100	2000	CHINS0129
		30G x 1/2"	12,7 mm	100	2000	CHINS0130
		30G x 5/16"	8 mm	100	2000	CHINS01308
		31G x 4/16"	6 mm	100	2000	CHINS01316
		Luer		100	2000	CHINS01
INSULIN U100	0,5 ml	29G x 1/2"	12,7 mm	100	2000	CHINS00529
		30G x 1/2"	12,7 mm	100	2000	CHINS00530
		31G x 5/16"	8 mm	100	2000	CHINS005318
		31G x 4/16"	6 mm	100	2000	CHINS005316
	0,3 ml	29G x 1/2"	12,7 mm	100	2000	CHINS00329
		30G x 1/2"	12,7 mm	100	2000	CHINS00330
		31G x 5/16"	8 mm	100	2000	CHINS003318
		31G x 4/16"	6 mm	100	2000	CHINS003316

CHIRANA® Tuberculin syringes

Type	Size	Size G/"	Size metric	Pieces per Box	Pieces per Case	Item Code
TUBERCULIN	1 ml	29G x 1/2"	12,7 mm	100	2000	CHTUB0129
		Luer		100	2000	CHTUB01

- ✓ possibility of 10 pcs POLY BAG packaging
- ✓ individual OEM manufacturing in reasonable quantity available



NEW



MEDOJECT fine® Insulin Pen needles

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ compatible with most of type A Insulin pens
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ shorter needles reduce risk of a painful injection into the muscle
- ✓ comfortable injection experience
- ✓ sterile, individually packed

Size	Color code	Pieces per Box	Pieces per Case	Item Code
29G x10 mm	yellow	100	2000	CHPN2910
29G x 12,7 mm	pink	100	2000	CHPN29127
30G x 6 mm	dark blue	100	2000	CHPN3006
30G x 8 mm	light blue	100	2000	CHPN3008
31G x 4 mm	green	100	2000	CHPN3104
31G x 5 mm	dark violet	100	2000	CHPN3105
31G x 6 mm	dark blue	100	2000	CHPN3106
31G x 8 mm	light blue	100	2000	CHPN3108
32G x 4 mm	green	100	2000	CHPN3204
32G x 6 mm	dark blue	100	2000	CHPN3206

- ✓ intended for Insulin application via most of type A Insulin pens
- ✓ made of polypropylene hub and stainless steel tube
- ✓ thin wall tube material





MEDOJECT® Hypodermic needles

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ precision lancet point grinding, silicone treated for smooth penetration
- ✓ color coded transparent hub for flow visual control
- ✓ included all metric and gauge measurements on all packages
- ✓ mainly thin wall tube material
- ✓ made of polypropylene hub and stainless steel tube
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ 100 pcs per box, 4000 pcs per case
- ✓ sterile, individually blister packed

Size metric	Size G/"	Color code	Item Code	Size metric	Size G/"	Color code	Item Code
0,30 x 12 mm	30 G x 1/2"	light yellow	CH30012	0,70 x 40 mm	22 G x 1 1/2"	black	CH22112
0,40 x 12 mm	27 G x 1/2"	grey	CH27012	0,70 x 50 mm	22 G x 2"	black	CH22200
0,40 x 20 mm	27 G x 4/5"	grey	CH27045	0,80 x 16 mm	21 G x 5/8"	green	CH21058
0,40 x 20 mm	27 G x 3/4"	grey	CH27034	0,80 x 25 mm	21 G x 1"	green	CH21100
0,40 x 40 mm	27 G x 1 1/2"	grey	CH27112	0,80 x 30 mm	21 G x 1 1/4"	green	CH21114
0,45 x 12 mm	26 G x 1/2"	brown	CH26012	0,80 x 40 mm	21 G x 1 1/2"	green	CH21112
0,45 x 16 mm	26 G x 5/8"	brown	CH26058	0,80 x 50 mm	21 G x 2"	green	CH21200
0,45 x 25 mm	26 G x 1"	brown	CH26100	0,90 x 25 mm	20 G x 1"	yellow	CH20100
0,50 x 16 mm	25 G x 5/8"	orange	CH25058	0,90 x 30 mm	20 G x 1 1/4"	yellow	CH20114
0,50 x 20 mm	25 G x 4/5"	orange	CH25045	0,90 x 40 mm	20 G x 1 1/2"	yellow	CH20112
0,50 x 25 mm	25 G x 1"	orange	CH25100	0,90 x 50 mm	20 G x 2"	yellow	CH20200
0,50 x 40 mm	25 G x 1 1/2"	orange	CH25112	1,10 x 25 mm	19 G x 1"	ivory	CH19100
0,55 x 25 mm	24 G x 1"	purple	CH24100	1,10 x 30 mm	19 G x 1 1/4"	ivory	CH19114
0,55 x 38 mm	24 G x 1 1/2"	purple	CH24112	1,10 x 40 mm	19 G x 1 1/2"	ivory	CH19112
0,60 x 16 mm	23 G x 5/8"	blue	CH23058	1,10 x 50 mm	19 G x 2"	ivory	CH19200
0,60 x 25 mm	23 G x 1"	blue	CH23100	1,20 x 25 mm	18 G x 1"	pink	CH18100
0,60 x 30 mm	23 G x 1 1/4"	blue	CH23114	1,20 x 40 mm	18 G x 1 1/2"	pink	CH18112
0,60 x 40 mm	23 G x 1 1/2"	blue	CH23112	1,20 x 50 mm	18 G x 2"	pink	CH18200
0,60 x 60 mm	23 G x 2 1/4"	blue	CH23214	1,50 x 50 mm	17 G x 2"	red purple	CH17200
0,70 x 25 mm	22 G x 1"	black	CH22100	1,60 x 25 mm	16 G x 1"	white	CH16100
0,70 x 30 mm	22 G x 1 1/4"	black	CH22114	1,60 x 40 mm	16 G x 1 1/2"	white	CH16112
0,70 x 35 mm	22 G x 1 2/5"	black	CH22125	1,80 x 40 mm	15 G x 1 1/2"	blue grey	CH15112

- ✓ intended for intramuscular, subcutaneous, intravenous and intra-arterial application
- ✓ possibility of short bevel lancet point grinding



CHIRAVAC™ Blood Collection needles

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ screw-on needle hub compatible with all available holders
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ easy-to-open packaging and manipulation
- ✓ comfortable collection experience
- ✓ 100 pcs per box, 1000 pcs per case
- ✓ sterile, individually packed

Size metric	Size G/"	Color code	Item Code	Size metric	Size G/"	Color code	Item Code
0,60 x 38 mm	23G x 1 1/2"	blue	CHBCN23112	0,60 x 25 mm	23G x 1"	blue	CHBCN23100
0,70 x 38 mm	22G x 1 1/2"	black	CHBCN22112	0,70 x 25 mm	22G x 1"	black	CHBCN22100
0,80 x 38 mm	21G x 1 1/2"	green	CHBCN21112	0,80 x 25 mm	21G x 1"	green	CHBCN21100
0,90 x 38 mm	20G x 1 1/2"	yellow	CHBCN20112	0,90 x 25 mm	20G x 1"	yellow	CHBCN20100
1,20 x 38 mm	18G x 1 1/2"	pink	CHBCN18112	1,20 x 25 mm	18G x 1"	pink	CHBCN18100

- ✓ intended for blood collection with vacuum tubes
- ✓ made of polypropylene hub and stainless steel tube
- ✓ thin wall tube material
- ✓ color coded needle hub



MEDOJECT® SAFETY Hypodermic needles

Size metric	Size G/"	Color code	Item Code	Size metric	Size G/"	Color code	Item Code
0,60 x 30 mm	23 G x 1 1/4"	blue	CHS23114	0,80 x 40 mm	21 G x 1 1/2"	green	CHS21112
0,70 x 40 mm	22 G x 1 1/2"	black	CHS22112	0,90 x 40 mm	20 G x 1 1/2"	yellow	CHS20112

- ✓ integrated safety feature for one-handed manipulation



MEDOJECT® Blunt fill needles

Size metric	Size G/''	Color code	Item Code
1,20 x 40 mm BLUNT FILL	18 G x 1 1/2''	pink	CH18112BF
1,20 x 40 mm BLUNT FILL w/ FILTER	18 G x 1 1/2''	pink	CH18112F
1,20 x 50 mm BLUNT FILL w/FILTER	18 G x 2''	pink	CH18200F

- ✓ 5µm FILTER for aspiration filtration procedure
- ✓ particles risk elimination during penetration



INOX® Ophthalmic needles

Size metric	Size G/''	Color code	Item Code STRAIGHT	Item Code BENT
0,40 x 20 mm	27 G x 4/5''	grey	CHOP27045S	CHOP27045B
0,60 x 20 mm	23 G x 4/5''	blue	CHOP23045S	CHOP23045B

- ✓ intended for eye irrigation
- ✓ end of needle free of burrs

MEDOJECT® Special Hypodermic needles

Size metric	Size G/''	Color code	Item Code
0,60 x 80 mm	23 G x 3 1/8''	blue	CH23018
0,80 x 100 mm	21 G x 4''	green	CH21400
0,80 x 120 mm	21 G x 4 3/4''	green	CH21434
0,90 x 70 mm	20 G x 2 4/5''	yellow	CH20245
1,20 x 70 mm	18 G x 2 4/5''	pink	CH18245
1,20 x 100 mm	18 G x 4''	pink	CH18400
1,60 x 70 mm	16 G x 2 4/5''	white	CH16245
1,80 x 70 mm	15 G x 2 4/5''	blue grey	CH15245
2,00 x 70 mm	14 G x 2 4/5''	pale green	CH14245
2,00 x 80 mm	14 G x 3 1/8''	pale green	CH14318
2,00 x 100 mm	14 G x 4''	pale green	CH14400
2,00 x 120 mm	14 G x 4 3/4''	pale green	CH14434

- ✓ intended for wide special application range



CHIRAFLEX® IV Cannulas | CHIRAFLEX® SAFETY IV Cannulas

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ intended for peripheral vascular application
- ✓ compatible with all standard medical devices for wide application
- ✓ catheter made from Teflon material, available in SAFETY version
- ✓ 50 pcs per box, 1000 pcs per case

Size	Diameter (mm)	Length (mm)	Flowrate (ml/min)	Color code	Item Code WITH PORT WITH WINGS	Item Code WITHOUT PORT WITH WINGS	Item Code WITHOUT PORT WITHOUT WINGS	Item Code SAFETY WITH PORT WITH WINGS
26 G	0,6	19	17	purple	CHFX0102619	CHFX0202619	CHFX0302619	CHFXS0102619
24 G	0,7	19	23	yellow	CHFX0102419	CHFX0202419	CHFX0302419	CHFXS0102419
22 G	0,9	25	36	blue	CHFX0102225	CHFX0202225	CHFX0302225	CHFXS0102225
20 G	1,1	32	56	pink	CHFX0102032	CHFX0202032	CHFX0302032	CHFXS0102032
18 G	1,3	45	90	green	CHFX0101845	CHFX0201845	CHFX0301845	CHFXS0101845
17 G	1,5	45	142	white	CHFX0101745	CHFX0201745	CHFX0301745	CHFXS0101745
16 G	1,7	45	200	grey	CHFX0101645	CHFX0201645	CHFX0301645	CHFXS0101645
14 G	2,1	45	305	orange	CHFX0101445	CHFX0201445	CHFX0301445	CHFXS0101445

CHIRAFLEX® Stylets

- ✓ compatible with IV Cannulas and SAFETY IV Cannulas
- ✓ color coded as per IV Cannula color coding

Size	Color Code	Pieces per Box	Pieces per Case	Item Code
22 G	blue	50	1000	CHFX20022
20 G	pink	50	1000	CHFX20020
18 G	green	50	1000	CHFX20018
17 G	white	50	1000	CHFX20017
16 G	grey	50	1000	CHFX20016
14 G	orange	50	1000	CHFX20014

- ✓ intended for longer interruptions in infusion
- ✓ screw onto needle adapter



CHIRAPLUS™ Infusion set

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ all sets includes 15 µm filter, ventilation port and luer lock adapter
- ✓ CHIRAPLUS P intended for pressure infusion, does not include flash ball
- ✓ CHIRAPLUS G intended for gravity infusion, includes flash ball
- ✓ all sets have 150 cm tubing, unless specified in the product item name

Item	Chamber Volume	Packaging type	Pieces per Case	Item Code
CHIRAPLUS P	10,5 ml	blister	100	CHP336
CHIRAPLUS P DEHP free	10,5 ml	blister	100	CHP536
CHIRAPLUS P 12	12,0 ml	blister	100	CHP306
CHIRAPLUS P 12 DEHP free	12,0 ml	blister	100	CHP506
CHIRAPLUS P 190	10,5 ml	blister	100	CHP339
CHIRAPLUS P 190 SAFETY	10,5 ml	blister	100	CHP338
CHIRAPLUS P DEHP Free Autofusion 180	15,0 ml	blister	300	CHP576
CHIRAPLUS G	10,5 ml	blister	100	CHG331
CHIRAPLUS G	10,5 ml	poly bag	500	CHG317
CHIRAPLUS G DEHP free	10,5 ml	blister	100	CHG531
CHIRAPLUS G 12	12,0 ml	blister	100	CHG301
CHIRAPLUS G 12 DEHP free	12,0 ml	blister	100	CHG501
CHIRAPLUS G	8,5 ml	poly bag	500	CHG318
CHIRAPLUS G DEHP free	8,5 ml	poly bag	500	CHG518
CHIRAPLUS G with needle	8,5 ml	poly bag	500	CHG319
CHIRAPLUS G with metal spike	8,5 ml	poly bag	500	CHG316

CHIRAHEM™ Transfusion set

Item	Chamber Volume	Packaging type	Pieces per Case	Item Code
CHIRAHEM	17,0 ml	blister	100	CHT314
CHIRAHEM DEHP Free	17,0 ml	blister	100	CHT514
CHIRAHEM	17,0 ml	poly bag	300	CHT320
CHIRAHEM with needle	17,0 ml	poly bag	300	CHT321
CHIRAHEM with metal spike	17,0 ml	poly bag	300	CHT322

- ✓ intended for blood transfusion administration
- ✓ all sets have 200 µm blood filter



CHIRAPLUS™ Spikes | CHIRAPLUS™ Injection ports | CHIRAPLUS™ Stoppers

- ✓ spikes intended for comfortable and safe drug manipulation
- ✓ ports offer additional drug administration when applied with other devices
- ✓ combistopper and stopper give a short-term option for storage or transfer of solution

Item	Pieces per Box	Pieces per Case	Item Code
Infusion spike RED air 0,2 µm + 5 µm solution filter	25	400	CHIS01
Infusion spike BLUE air 1,2 µm + 5 µm solution filter	25	400	CHIS02
Infusion spike GREEN air 1,2 µm filter	25	400	CHIS03
Transflo spike	100	1000	CHTS01
Needle free injection port	100	2400	CHNF01
Injection port	100	5000	CHIP01
Combistopper RED	100	5000	CHCS01
Combistopper BLUE	100	5000	CHCS02
Combistopper WHITE	100	5000	CHCS03
Stopper Luer lock	100	5000	CHLL01

- ✓ intended for range of application in IV administration



CHIRAFLEX® Scalp vein sets

Size	Diameter	Color code	Pieces per Box	Pieces per Case	Item Code
27 G	0,4 mm	grey	100	2000	CH260001127
25 G	0,5 mm	orange	100	2000	CH260001125
23 G	0,6 mm	blue	100	2000	CH260001123
22 G	0,7 mm	black	100	2000	CH260001122
21 G	0,8 mm	green	100	2000	CH260001121
20 G	0,9 mm	yellow	100	2000	CH260001120
19 G	1,1 mm	ivory	100	2000	CH260001119
18 G	1,2 mm	pink	100	2000	CH260001118

- ✓ stainless steel needle with length of 20 mm (4/5"), 30 cm plastic tubing
- ✓ intended for shorter time infusion or pediatric use



CHIRALINE™ Extension lines

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ available extension tubing in 45 cm or 180 cm
- ✓ luer lock ends available
- ✓ compatible with any available medical device

Item	Pieces per Box	Pieces per Case	Item Code
CHIRALINE 450	120	960	CHL318045
CHIRALINE 1800	85	680	CHL318180

- ✓ intended for extension of infusion or transfusion range of application, OD 3 mm x ID 1,95 mm



CHIRAWAY™ 3-way stopcock

- ✓ available in standard or lipid resistant material
- ✓ color coding differentiation
- ✓ compatible with any available medical device

Item	Pieces per Box	Pieces per Case	Item Code
CHIRAWAY Standard	50	1000	CHW001
CHIRAWAY Lipid resistant	50	1000	CHW002

- ✓ intended for infusion or transfusion range of application



CHIRASKIN™ Gloves

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ non-sterile, intended for examination
- ✓ high elasticity, high strength, excellent grip sensitivity
- ✓ safe protection, good touch sense

Size	Pieces per Box	Pieces per Carton	Item Code LATEX POWDERED	Item Code LATEX POWDER FREE	Item Code NITRIL POWDER FREE
Gloves Small	100	1000	CHL01S	CHLF01S	CHNF01S
Gloves Medium	100	1000	CHL01M	CHLF01M	CHNF01M
Gloves Large	100	1000	CHL01L	CHLF01L	CHNF01L

Examination products

Item	Description	Pieces per Box	Pieces per Carton	Item Code
Face mask TIE ON	3-ply, non-woven	50	1000	CHFM01
Face mask EAR LOOP	3-ply, non-woven	50	1000	CHFM02
Bouffant cap 21"	non-woven	100	1000	CHC01
Bouffant cap 24"	non-woven	100	1000	CHC02
Shoecovers	CPE, blue	100	2000	CHSC01
Bedcover	CPE, blue	10	100	CHBC01

- ✓ intended for variety of application range
- ✓ support for safe and clean healthcare

CTI - MD|ENG|03|2018

MEDICAL DEVICES

 **CHIRANA T. Injecta**

CHIRANA T. Injecta, a.s., Nám. Dr. Schweitzera 194
916 01 Stará Turá, Slovak Republic

www.t-injecta.sk