

MEDOJECT® Hypodermic needles



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ precision lancet point grinding, silicone treated for smooth penetration
- ✓ colored coded transparent hub for flow visual control
- ✓ included all metric and gauge measurements on all packages
- ✓ mainly thin wall tube material
- ✓ made of polypropylene hub and stainless steel tube
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ 100 pcs per box, 4000 pcs per case sterile, individually blister packed

Size metric	Size G/''	Color code	Item Code
0,30 x 12 mm	30 G x 1/2''	light yellow	CH30012
0,40 x 12 mm	27 G x 1/2''	grey	CH27012
0,40 x 20 mm	27 G x 4/5''	grey	CH27045
0,40 x 20 mm	27 G x 3/4''	grey	CH27034
0,40 x 40 mm	27 G x 1 1/2''	grey	CH27112

MEDOJECT® Hypodermic needles

Size metric	Size G/''	Color code	Item Code
0,45 x 12 mm	26 G x 1/2''	brown	CH26012
0,45 x 16 mm	26 G x 5/8''	brown	CH26058
0,45 x 25 mm	26 G x 1''	brown	CH26100
0,50 x 16 mm	25 G x 5/8''	orange	CH25058
0,50 x 20 mm	25 G x 4/5''	orange	CH25045
0,50 x 25 mm	25 G x 1''	orange	CH25100
0,50 x 40 mm	25 G x 1 1/2''	orange	CH25112
0,55 x 25 mm	24 G x 1''	purple	CH24100
0,55 x 38 mm	24 G x 1 1/2''	purple	CH24112
0,60 x 16 mm	23 G x 5/8''	blue	CH23058
0,60 x 25 mm	23 G x 1''	blue	CH23100
0,60 x 30 mm	23 G x 1 1/4''	blue	CH23114
0,60 x 40 mm	23 G x 1 1/2''	blue	CH23112
0,60 x 60 mm	23 G x 2 1/4''	blue	CH23214
0,70 x 25 mm	22 G x 1''	black	CH22100
0,70 x 30 mm	22 G x 1 1/4''	black	CH22114
0,70 x 35 mm	22 G x 1 2/5''	black	CH22125
0,70 x 40 mm	22 G x 1 1/2''	black	CH22112
0,70 x 50 mm	22 G x 2''	black	CH22200
0,80 x 16 mm	21 G x 5/8''	green	CH21058
0,80 x 25 mm	21 G x 1''	green	CH21100

MEDOJECT® Hypodermic needles

Size metric	Size G/''	Color code	Item Code
0,80 x 30 mm	21 G x 1 1/4''	green	CH21114
0,80 x 40 mm	21 G x 1 1/2''	green	CH21112
0,80 x 50 mm	21 G x 2''	green	CH21200
0,90 x 25 mm	20 G x 1''	yellow	CH20100
0,90 x 30 mm	20 G x 1 1/4''	yellow	CH20114
0,90 x 40 mm	20 G x 1 1/2''	yellow	CH20112
0,90 x 50 mm	20 G x 2''	yellow	CH20200
1,10 x 25 mm	19 G x 1''	ivory	CH19100
1,10 x 30 mm	19 G x 1 1/4''	ivory	CH19114
1,10 x 40 mm	19 G x 1 1/2''	ivory	CH19112
1,10 x 50 mm	19 G x 2''	ivory	CH19200
1,20 x 25 mm	18 G x 1''	pink	CH18100
1,20 x 40 mm	18 G x 1 1/2''	pink	CH18112
1,20 x 50 mm	18 G x 2''	pink	CH18200
1,50 x 50 mm	17 G x 2''	red purple	CH17200
1,60 x 25 mm	16 G x 1''	white	CH16100
1,60 x 40 mm	16 G x 1 1/2''	white	CH16112
1,80 x 40 mm	15 G x 1 1/2''	blue grey	CH15112

✓ intended for intramuscular, subcutaneous, intravenous and intra-arterial application

✓ possibility of short bevel lancet point grinding

