

MEDOJECT® Special Hypodermic needles



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ intended for wide special application range
- ✓ possibility of short bevel lancet point grinding
- ✓ availability of blunt fill or blunt fill with filter needle
- ✓ other not listed sizes available upon request

Size metric	Size G/''	Color code	Item Code
0,60 x 80 mm	23 G x 3 1/8''	blue	CH23318
0,80 x 100 mm	21 G x 4''	green	CH21400
0,80 x 120 mm	21 G x 4 3/4''	green	CH21434
0,90 x 70 mm	20 G x 2 4/5''	yellow	CH20245
1,20 x 70 mm	18 G x 2 4/5''	pink	CH18245
1,20 x 100 mm	18 G x 4''	pink	CH18400

MEDOJECT® Special Hypodermic needles

Size metric	Size G/''	Color code	Item Code
1,60 x 70 mm	16 G x 2 4/5''	white	CH16245
1,80 x 70 mm	15 G x 2 4/5''	blue grey	CH15245
2,00 x 70 mm	14 G x 2 4/5''	pale green	CH14245
2,00 x 80 mm	14 G x 3 1/8''	pale green	CH14318
2,00 x 100 mm	14 G x 4''	pale green	CH14400
2,00 x 120 mm	14 G x 4 3/4''	pale green	CH14434