

MEDOJECT® fine Insulin Pen needles



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ compatible with most of type A Insulin pens
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ shorter needles reduce risk of a painful injection into the muscle
- ✓ comfortable injection experience
- ✓ sterile, individually packed

| Size | Color code | Pieces per Box | Pieces per Case | Item Code |
|---------------|-------------|----------------|-----------------|-----------|
| 29G x10 mm | yellow | 100 | 2000 | CHPN2910 |
| 29G x 12,7 mm | pink | 100 | 2000 | CHPN29127 |
| 30G x 6 mm | dark blue | 100 | 2000 | CHPN3006 |
| 30G x 8 mm | light blue | 100 | 2000 | CHPN3008 |
| 31G x 4 mm | green | 100 | 2000 | CHPN3104 |
| 31G x 5 mm | dark violet | 100 | 2000 | CHPN3105 |

MEDOJECT® fine Insulin Pen needles

| Size | Color code | Pieces per Box | Pieces per Case | Item Code |
|------------|-------------|----------------|-----------------|-----------|
| 31G x 6 mm | dark blue | 100 | 2000 | CHPN3106 |
| 31G x 8 mm | light blue | 100 | 2000 | CHPN3108 |
| 32G x 4 mm | green | 100 | 2000 | CHPN3204 |
| 32G x 5 mm | dark violet | 100 | 2000 | CHPN3205 |
| 32G x 6 mm | dark blue | 100 | 2000 | CHPN3206 |

- ✓ intended for Insulin application via most of type A Insulin pens
- ✓ made of polypropylene hub and stainless steel tube
- ✓ thin wall tube material